

## PRODUCT DATA

# TSPFL

<u>SUBSTANCE</u>		<u>THICKNESS</u>	
Adhesive	35gsm		
Tissue	17gsm	Adhesive & Carrier	85 micron
Adhesive	35gsm		
Liner	65gsm	Liner	55 micron

### PRODUCT DESCRIPTION

A permanent tape with a tissue carrier. The adhesive is an aggressive rubber based system and the liner a two side release system which overlaps at each edge to give a 'finger-lift' for easy liner removal.

### APPLICATIONS

The original and industry standard finger-lift tissue tape, DST/F has been used for continuous stationery, mailing envelopes and padded envelopes, studio/commercial artwork, photographic, Point of Sale and splicing. Especially useful for many products produced by print finishers.

### ADHESIVE PROPERTIES

		<u>20 min</u>	<u>24 hour</u>	<u>TEST METHOD</u>
Adhesion to steel	open side	1300 N/M	1300 N/M	BS7116:1990 - D
	closed side	1300 N/M	1300 N/M	BS7116:1990 - D
Shear Adhesion	open side	100 hours		BS7116:1990 - B
	closed side	100 hours		BS7116:1990 - B
Release of liner from adhesive		25 N/M		BS7116:1990 - E

### THERMAL PROPERTIES

Recommended service temperature range	0°C to 70°C
Application temperature range	10°C to 50°C

### SHELF LIFE

Minimum shelf life when stored in closed cartons away from direct sources of heat and sunlight	2 Years
--	---------

### HEALTH AND SAFETY INFORMATION

Adhesive conforms to F&DA Regulation 175-105